

BULLSEYE GRID



- Targets injuries, trigger or acu points for conditions.
- Reduces pain.
- Reduces numbness.
- Finds blockages (impedance).
- Shows progress in readings.

BILATERAL SPINAL



- Targets Bladder Meridian.
- Benefits all conditions.
- Regulates central nervous system.
- Promotes pain relief.
- Improves sleep, energy & immune function.
- Enhances emotional balance.
- Improves urinary health.

PLATES



- Targets any condition.
- Localizes directly to point of pain or problem with:
 - Origin/insertion.
 - Myofascial trigger points.
 - Acupressure points.
 - Neurolymph reflex points.
 - Dermatomes, etc.

AURICULAR



- Targets any condition & whole body.
- Ear/hearing issues.
- Shoulder issues.
- Wellbeing.
- Pain relief.
- Addiction issues.
- For specific conditions use a plate on the problem area and stimulate targeted points on the ear.

BILATERAL HANDS & FEET (IMMERSION)



- Provides charge and support for the whole body.
- Mild detox.
- Good systemic for chronic issues.
- Arthritis.
- Wellbeing.

EAR CLIPS



- Targets brain, head, neck and memory issues.
- Promotes relaxation.
- Stands alone or can be combined with other protocols.
- Ensure that patient/client is comfortable with cranial stimulation.
- Contraindicated in patients/clients with a history of seizures.
- Provide a blanket for warmth.

SP 6 & ST 36



- SP6 - for urogenital, muscles, legs, edema, immune system, insomnia.
- ST36 - for gastrointestinal, ulcers, reflux, constipation, enhancing energy.
- SP6 & ST36 - for gastrointestinal & urogenital issues.

Y-PROBE MASSAGE



- Targets inflammation.
- Used for dispersion of cellular waste after Bullseye Grid or Plates.
- Increases range of motion (flexibility).
- Promotes quick pain relief.

SEALING



- Targets the protective layer of electricity on the skin.
- Replaces skin static electricity after Bullseye Grid, Spinal Protocols, etc.

HEADBAND



- Targets brain, head, neck and emotional-related issues.
- Stands alone or can be combined with other protocols.
- Ensure that patient/client is comfortable with cranial stimulation.
- Contraindicated in patient/client with a history of seizures.
- Provide a blanket for warmth.

GOVERNING VESSEL



- Targets the Master.
- Benefits all conditions.
- Improves spine & back health.
- Contributes to stress reduction & balancing energies.
- Improves mental clarity & focus.
- Improves circulation.
- Enhances sexual health.

GUIDED MEDITATION



- Enhances wellbeing & physical health.
- Promotes better sleep.
- Helps with pain management.
- Boosts self-awareness & creativity.
- Stands alone or can be combined with other protocols.
- Provide a blanket for warmth.
- Continued healing with home practice.

ODONTON



- Targets any condition & whole body.
- Jaw or teeth-related issues.
- Wellbeing.
- Gastrointestinal track issues.
- Limbic system conditions.
- Sensation usually reduces by the third stimulation on the reactive tooth root.
- For specific conditions use a plate on the problem area and stimulate targeted tooth-organ access point.

LI 4 & LI 11



- LI4 - for head, neck, face, mouth, nasal issues.
- LI11 - for arm pain, fever, inflammation, allergies, skin conditions, constipation.
- LI4 & LI4 - clears head pathogens, induces sweating, opens channels.
- Targets large intestine issues.
- Promotes whole body pain relief.

REFLEXOLOGY



- Targets any condition & whole body.
- Hand- or foot-related issues.
- Wellbeing.
- Pain relief.
- Improves circulation, detoxification.
- For specific conditions use a plate on the problem area and stimulate targeted hand/foot reflexology access point.

LYMPH DRAINAGE



- Enhances circulation.
- Improves tissue environment.
- Supports whole body detoxification.
- Improves white blood cells and immune function.

ENDORPHIN PUMP



- Targets adrenal and vagus nerve points.
- Stimulates mood and clarity.
- Creates calm.
- May cause vertigo feeling during stimulation.

BEAUTY BASIC



- Softer smoother skin.
- Lifting and tightening benefits.
- Cumulative effect.